



**TRACK
YOUR
VALUABLE
ASSETS**

**With QTRMaster
RFID Solution**

Loan Filtering

The Basics

2024-02-01


Loan Filtering

Release 1.4.0-68:

In order to evaluate the performance of “store” room locations as well as quickly reference all transactions that are recorded by borrowers using the ‘touch’ terminals within your equipment rooms, we have augmented our software to allow for filters within the main “LOAN” page: These additional filters are “Location In”, “Location Out” & “Last Seen”

Location ‘IN’

Selecting the ‘Location In’ filter will also list two extra fields “Item Damaged” and “Item Returned”



These options relate to the touch screen terminal usage at the equipment rooms and will display all transactions where the borrowers engaged the touch screens to report items which are not ‘seen by the reader’ as returned with their original transaction, however HAVE been returned.

Loan Filtering

Location In: Choose the 'Store' (room) reader you wish to monitor and the QTRMaster system will display all activities 'direction in', past and present, from selected 'Store' (room) reader.

Loans
IDALLOGY

- Dashboard
- People
- Items
- Loans
- Tags
- Locations
- Maintenance Needed
- Audit
- Configs
- Item Types
- Roles
- Rules
- Maintenance Types
- Reason Choices
- Settings
- Administrators
- Scanners
- Imports
- Provisions
- Brandings
- Licensings
- Services

Location In
Lab
ADD FILTER
EXPORT

Borrower	Item	Date Loaned ↓	Date Returned	Return expected time	Returned by	Location In	Last Seen Location	Date Seen	Direction Seen	Violation
Smith,David	AED yellow (36) [245]	2024-01-10, 7:52:28 p.m.	2024-01-10, 7:52:40 p.m.		Smith,David	Lab				
Smith,David	C8 () [2]	2024-01-10, 7:52:28 p.m.	2024-01-10, 7:52:40 p.m.		Smith,David	Lab				Not Allowed
Nobody	AED yellow (36) [245]	2023-10-27, 7:38:08 a.m.	2023-10-27, 7:38:15 a.m.		Nobody	Lab				
Smith,David	C8 () [2]	2023-10-14, 5:28:45 p.m.	2023-10-14, 5:29:14 p.m.		Smith,David	Lab				Not Allowed
Smith,David	AED yellow (36) [245]	2023-10-14, 5:28:45 p.m.	2023-10-14, 5:29:14 p.m.		Smith,David	Lab				
Nobody	Radio (2nd)	2023-10-07, 10:07:12 a.m.	2024-01-10, 8:07:58 p.m.		Nobody	Lab				Late Return
Smith,David	AED (1234)	2023-08-18, 8:18:58 p.m.	2023-08-18, 8:19:04 p.m.		Smith,David	Lab				
Smith,David	Cylinder (223)	2023-08-18, 8:18:58 p.m.	2023-08-18, 8:19:04 p.m.		Smith,David	Lab				Wrong Location
Smith,David	Ear Protection 2 ()	2023-08-18, 8:18:58 p.m.	2023-08-18, 8:19:04 p.m.		Smith,David	Lab				Late Return
Nobody	Holster Pistol (44111)	2023-08-18, 8:07:25 p.m.	2023-08-18, 8:18:58 p.m.		Nobody	Lab				Late Return

Rows per page: 10 ▾ | 1-10 of 8905 | 1 2 ... 891 NEXT >

Loan Filtering

Location Out: Choose the 'Store' (room) reader you wish to monitor and the QTRMaster system will display all activities 'direction out', past and present, from selected 'Store' (room) reader

The screenshot displays the IDALLOGY Loans management interface. The left sidebar contains navigation options: Dashboard, People, Items, Loans (selected), Tags, Locations, Maintenance Needed, Audit, Configs (Item Types, Roles, Rules, Maintenance Types, Reason Choices), Settings (Administrators, Scanners, Imports, Provisions, Brandings, Licensings, Services). The main content area shows a table of loan records with a filter for 'Location Out'. The table has columns: Borrower, Item, Date Loaned, Date Returned, Return expected time, Returned by, Location In, Last Seen Location, Date Seen, Direction Seen, and Violation. The table contains 12 rows of data, with violations such as 'Not Allowed', 'Not out', and 'Late Return' indicated by colored buttons. The bottom right of the table shows pagination: 'Rows per page: 10', '1-10 of 17802', and 'NEXT >'.

Borrower	Item	Date Loaned ↓	Date Returned	Return expected time	Returned by	Location In	Last Seen Location	Date Seen	Direction Seen	Violation
Smith,David	AED yellow (36) [245]	2024-01-10, 7:52:28 p.m.	2024-01-10, 7:52:40 p.m.		Smith,David	Lab				
Smith,David	C8 () [2]	2024-01-10, 7:52:28 p.m.	2024-01-10, 7:52:40 p.m.		Smith,David	Lab				Not Allowed
Smith,David	C8 () [2]	2024-01-10, 7:52:09 p.m.	2024-01-10, 7:52:09 p.m.							Not out
Smith,David	AED yellow (36) [245]	2024-01-10, 7:52:09 p.m.	2024-01-10, 7:52:09 p.m.							Not out
Nobody	AED yellow (36) [245]	2023-10-27, 7:38:08 a.m.	2023-10-27, 7:38:15 a.m.		Nobody	Lab				
Smith,David	C8 () [2]	2023-10-14, 5:28:45 p.m.	2023-10-14, 5:29:14 p.m.		Smith,David	Lab				Not Allowed
Smith,David	AED yellow (36) [245]	2023-10-14, 5:28:45 p.m.	2023-10-14, 5:29:14 p.m.		Smith,David	Lab				
Nobody	Radio (2nd)	2023-10-07, 10:07:12 a.m.	2024-01-10, 8:07:58 p.m.		Nobody	Lab				Late Return
Nobody	AED yellow (36) [245]	2023-09-06, 11:04:58 a.m.	2023-09-06, 11:04:58 a.m.							Not out
Nobody	C8 () [2]	2023-09-06, 11:04:58 a.m.	2023-09-06, 11:04:58 a.m.							Not out

Loan Filtering

The screenshot shows the IDALLOGY Loans management interface. On the left is a navigation sidebar with categories like Dashboard, People, Items, Loans (selected), Tags, Locations, Maintenance Needed, Audit, Configs, Item Types, Roles, Rules, Maintenance Types, Reason Choices, Settings, Administrators, Scanners, Imports, Provisions, Brandings, Licensings, and Services. The main header includes the IDALLOGY logo and user profile icons. The main content area features a table of loan records with columns: Borrower, Item, Date Loaned, Date Returned, Return expected time, Returned by, Location In, Last Seen Location, Date Seen, Direction Seen, and Violation. A dropdown menu is open over the 'Borrower' column, showing 'Exit1' and 'Logistics' options. The table contains 10 rows of data with various violation statuses like 'Not Allowed', 'Not out', and 'Late Return'. At the bottom right, there is a pagination control showing 'Rows per page: 10', '1-10 of 17802', and page numbers '1', '2', '...', '1781', 'NEXT >'.

Borrower	Item	Date Loaned	Date Returned	Return expected time	Returned by	Location In	Last Seen Location	Date Seen	Direction Seen	Violation
Smith,David	AED yellow (36) [245]	2024-01-10, 7:52:28 p.m.	2024-01-10, 7:52:40 p.m.		Smith,David	Lab				
Smith,David	C8 () [2]	2024-01-10, 7:52:28 p.m.	2024-01-10, 7:52:40 p.m.		Smith,David	Lab				Not Allowed
Smith,David	C8 () [2]	2024-01-10, 7:52:09 p.m.	2024-01-10, 7:52:09 p.m.							Not out
Smith,David	AED yellow (36) [245]	2024-01-10, 7:52:09 p.m.	2024-01-10, 7:52:09 p.m.							Not out
Nobody	AED yellow (36) [245]	2023-10-27, 7:38:08 a.m.	2023-10-27, 7:38:15 a.m.		Nobody	Lab				
Smith,David	C8 () [2]	2023-10-14, 5:28:45 p.m.	2023-10-14, 5:29:14 p.m.		Smith,David	Lab				Not Allowed
Smith,David	AED yellow (36) [245]	2023-10-14, 5:28:45 p.m.	2023-10-14, 5:29:14 p.m.		Smith,David	Lab				
Nobody	Radio (2nd)	2023-10-07, 10:07:12 a.m.	2024-01-10, 8:07:58 p.m.		Nobody	Lab				Late Return
Nobody	AED yellow (36) [245]	2023-09-06, 11:04:58 a.m.	2023-09-06, 11:04:58 a.m.							Not out
Nobody	C8 () [2]	2023-09-06, 11:04:58 a.m.	2023-09-06, 11:04:58 a.m.							Not out

Choose the option from your 'Monitoring' readers (Exit Readers) and the system will show the last monitoring reader to have tracked your asset.